

## Indications of the Freka® Pexact, ENFIT



Long-term enteral nutrition especially:

- Oesophageal passage dysfunctions
  - Stenosis
  - Varices
  - Diverticulum
  - Corrosion burns
  - Radio- or chemotherapy
- To avoid ENT-tumor cell seeding
- After head- and mouth surgery
  - Only transnasal endoscopy
- Advantage:
  - Non-endoscopic removable
  - Improved durability of the Feeding Tube

## Ordering information

	Product	Packaging	Art.No.
Freka® Pexact II FR 15, ENFit Direct Puncture Gastrostomy		1 x 1	7601365
Gastropexy Device II		1 x 1	7601363
Gastropexy Device		1 x 1	7601364
Freka® GastroTube FR 15 Percutaneous Replacement System		1 x 1	7755648
Freka® Seldinger Wire Teflon-coated wire, 80 cm with flexible tube tip of 10 cm		1 x 1	7989911
<b>Accessories</b>			
		Packaging	Art.No.
Dressing Set acc. to Hermann		1 x 1	7901081
Erlanger Dressing Set		1 x 1	7981941
CD ROM Freka® Pexact Video about placement method		1 x 1	7317181

## Freka® Pexact II FR 15, ENFit

Direct Puncture Gastrostomy  
with safe and reliable Gastropexy



# Freka® Pexact, ENFit


Direct Puncture Gastrostomy with safe and reliable Gastropexy according to

Dr. Toshio Kadota & Dr. Fumiaki Ueno (Direct Puncture Method),  
 Dr. Masaki Funada (Gastropexy Device = Loop Fixture & Loop Fixture II), Japan and  
 Dr. Arno Dormann, Germany



## Gastropexy with innovative suturing method

- Improved fixation device enabling single-hand operation
- Operation of the device while monitoring the endoscopic image
- Safe fixation of stomach to abdominal wall
- Extra corporal suture technique
- No risk of peritoneal infection by direct access
- Short, straight stoma tracts
- Needle protection to avoid injuries



Insertion of Gastropexy Device (Loop fixture)



Overtaking thread with loop wire



Fixing thread




Gastropexy sutures

## Placement method

Freka® Pexact\*

with endoscopic control



\*For details of the procedure, refer to the instructions for use.

# Freka<sup>®</sup> Pexact, ENFit

Direct Puncture Gastrostomy with safe and reliable Gastropexy according to

Dr. Toshio Kadota & Dr. Fumiaki Ueno (Direct Puncture Method),  
 Dr. Masaki Funada (Gastropexy Device = Loop Fixture & Loop Fixture II), Japan and  
 Dr. Arno Dormann, Germany



## Gastropexy with innovative suturing method

- Improved fixation device enabling single-hand operation
- Operation of the device while monitoring the endoscopic image
- Safe fixation of stomach to abdominal wall
- Extra corporal suture technique
- No risk of peritoneal infection by direct access
- Short, straight stoma tracts
- Needle protection to avoid injuries



Insertion of Gastropexy Device (Loop fixture)




Overtaking thread with loop wire



Fixing thread



Gastropexy sutures



## Indications of the Freka<sup>®</sup> Pexact, ENFIT



Long-term enteral nutrition especially:

- Oesophageal passage dysfunctions
  - Stenosis
  - Varices
  - Diverticulum
  - Corrosion burns
  - Radio- or chemotherapy
- To avoid ENT-tumor cell seeding
- After head- and mouth surgery
  - Only transnasal endoscopy
- Advantage:
  - Non-endoscopic removable
  - Improved durability of the Feeding Tube