

# Ordering information

	Size in ml/mm	Packaging	Art.No.
<b>Enteral Dispensing Systems, ENFit</b>			
Freka® Connect, ENFit/ProNeo	10 ml	10 x 1	9000784XS
Freka® Connect, ENFit/ProNeo	10 ml	30 x 1	9000784S
Freka® Connect, ENFit/ProNeo	10 ml	100 x 1	9000784
Freka® Connect, ENFit/ProNeo	20 ml	10 x 1	9000785XS
Freka® Connect, ENFit/ProNeo	20 ml	30 x 1	9000785S
Freka® Connect, ENFit/ProNeo	20 ml	100 x 1	9000785
Freka® Connect, ENFit/ProNeo	60 ml	10 x 1	9000786XS
Freka® Connect, ENFit/ProNeo	60 ml	30 x 1	9000786S
Freka® Connect, ENFit/ProNeo	60 ml	50 x 1	9000786
Freka® Connect, ENFit	100 ml	10 x 1	9000787XS
Freka® Connect, ENFit	100 ml	30 x 1	9000787
Freka® Connect, ENFit/ProNeo, Eccentric	60 ml	50 x 1	9000789

## Low Dose Tip Syringes

Freka® Connect ENFit/ProNeo LDT	1 ml	30 x 1	9000788LDT
Freka® Connect ENFit/ProNeo LDT	2.5 ml	30 x 1	9000782LDT
Freka® Connect ENFit/ProNeo LDT	5 ml	30 x 1	9000783LDT

## Accessories

Syringe Adapter		100 x 1	7752014
Freka® Cap ENFit/ProNeo, sterile		30 x 1	7981402
Freka® Cap ENFit/ProNeo		1 x 30	7981403
Syringe Holder		1 x 10	7752020
ProNeo Multifunction Syringe Cap, (20 ml and 60 ml)		30 x 1	7752002
ProNeo Sampling Straw, sterile	50 mm	30 x 1	7752019
ProNeo Sampling Straw, sterile	100 mm	30 x 1	7752017
ProNeo Sampling Straw, sterile	150 mm	30 x 1	7752016
ProNeo Sampling Straw, non-sterile	100 mm	1 x 30	ENFT1001
ProNeo Sampling Straw, non-sterile	150 mm	1 x 30	7752013
ProNeo Sampling Cannula, sterile	50 mm	30 x 1	7752012
ProNeo MilkBottle Adapter		30 x 1	7752003
ProNeo Medication Adapter		30 x 1	7752004
VarioBottle Adapter		30 x 1	7752015
ProNeo Bottle Insert, non-removable	17-20,5 mm	30 x 1	9000793
ProNeo Bottle Insert, non-removable	12-15 mm	30 x 1	9000794
ProNeo Bottle Insert, non-removable	14-16,5 mm	30 x 1	9000795
ProNeo Bottle Insert, non-removable	19-23 mm	30 x 1	9000796
ProNeo Bottle Insert, non-removable	24-26 mm	30 x 1	9000797

The Freka® Connect Syringes are DEHP- and Latex-free.  
For more information on the ENFit syringes please contact your local Fresenius Kabi representative.  
ENFit® is a trademark of GEDSA.  
For more information visit the website [www.stayconnected.org](http://www.stayconnected.org).

# Freka® Connect, ENFit

## Dedicated Syringes

### for all Enteral Nutrition Applications



© Fresenius Kabi Deutschland GmbH. Reproduction and distribution - only with prior written authorization

7392891Z (12.22/AC)

# Freka® Connect, ENFit

Dedicated Syringes for all Enteral Nutrition Applications



## ENFit Centric Tip Syringes

The standard choice for enteral nutrition

- good visible scale
- high quality application haptic
- sterile, single packed
- various sizes: 10, 20, 60 and 100ml

## ENFit Eccentric Tip Syringes

Preferred solution for breast milk delivery:

- the eccentric tip improves the fat delivery by pump application compared to centric tips
- available in 20 and 60ml size

## ENFit Low Dose Tip Syringes

Precise application of small volumes

- good visible scale
- reduced dead space
- available in 1, 2.5 and 5ml size

**ENFit®**

**ENFit for safety** - All syringes are equipped with the safe enteral connector ENFit, according to the international safety standard ISO 80369-3. The ENFit system reduces the risk of an unintended misconnection. Therefore these syringes offer the highest safety in the clinical routine for patients and caregivers.

**ENFit = one enteral connector for all applications in enteral nutrition**

## Useful Syringe Accessories

Lower Contamination Risks by Using the right Accessories



### Adapter

- Easy and tight connection of two ENFit syringes to mix or transfer liquids



### Multifunction Syringe Cap

- With air filter to protect liquids and nutrition during gravity feeding
- Flexible in use, as it fits to 20 ml and 60 ml syringes



### Syringe Holder

- Syringe Holder for gravity feeding
- Parallel fixation of two ENFit syringes
- Flexible elastic strap for easy fixation



### Sampling Straw

- Clean and hygienic draw of liquids
- Protective cover to prevent kinking before usage
- Available in 5, 10 and 15 cm length



### Sampling Cannula

- Hygienic and precise draw of even small portions of liquids from glass ampoules or vials
- Integrated particle filter for additional safety
- Protective cover



### MilkBottle Adapter

- Ventilated adapter to draw milk or liquids out of a milk bottle
- Direct connection of an ENFit syringe



### Medication Adapter

- Vented, universal spike to draw medicine out of a sterile vial



### VarioBottle Adapter

- Adapter to draw liquids from a bottle or vial
- Flexible in use, as it fits to various bottle neck sizes



### Bottle Insert, non-removable

- For small bottles, e.g. medicine bottles
- To draw medicine out of a small bottle or vial
- Available in 5 different diameters



### Syringe Caps for ENFit syringes

- For a hygienic intermediate storage