Freka® PEG Pro with ENFit

... for patients without gastric motility disorders & without risk of aspiration.

Getting a PEG does not mean the end of a self-determined life. Rather, it is the opportunity to decisively improve the quality of life through a good nutritional status.

To meet the individual needs and find the most comfortable solution for the patients:

- Available in FR 12, 14, 16 and 20*
- Versatile usage
- Extension with an intestinal tube is possible
- Changeover to a balloon system, like Freka Belly Button

Several placement videos and a comprehensive digital aftercare training, can be explored under **www.enteral-nutrition.com**.

^{*} Freka PEG FR 9 and FR 15 are further available



Ordering information

	FR size	Packaging	Colour Code	Art.No.
Freka® PEG Pro, ENFit	FR 12	1 x 1		M90800346
	FR 14	1 x 1		M90800347
	FR 16	1 x 1		M90800349
	FR 20	1 x 1		M90800348
Freka® PEG, ENFit	FR 9	1 x 1		7755642
	FR 15	1 x 1		7755643
Freka® PEG gastric / Intestinal (PEG/J or JET-PEG), ENFit	FR 9/15	1 x 1		M90800351
	FR 9/16	1 x 1		M90800373
	FR 12/20	1 x 1		M90800350
			•	

	FR size	Packaging	Art.No.
Accessories			
Freka® Repair Sets for Freka® PEG Pro, ENFit	FR 12	1 x 1	M90800400
	FR 14	1 x 1	M90800401
	FR 16	1 x 1	M90800402
	FR 20	1 x 1	M90800403
Freka® Repair Sets for Freka® PEG, ENFit	FR 9	1 x 1	7981382
	FR 15	1 x 1	7981383
Dressing Set acc. to Hermann		1 x 1	7901081
Erlanger Dressing Set		1 x 1	7981941

	Size in ml	Packaging	Art.No.
Enteral Dispensing Systems, ENFit			
Freka Connect ENFit/ProNeo	20 ml	100 x 1	9000785
	20 ml	30 x 1	9000785S
	20 ml	10 x 1	9000785XS
	60 ml	50 x 1	9000786
	60 ml	30 x 1	9000786S
	60 ml	10 x 1	9000786XS
	100 ml	30 x 1	9000787
	100 ml	10 x 1	9000787XS



Fresenius Kabi AG Else-Kröner-Straße 1 61352 Bad Homburg v.d.H. Germany Phone: +49 (0) 6172 686-0 www.fresenius-kabi.com





Freka® PEG Pro, ENFit

With safety comes comfort



Freka® PEG Pro

Engineered for safe use at every single step

Percutaneous Endoscopic Gastrostomy* The choice for long-term enteral nutrition (> 4 weeks)

Safe usage even when pressed for time









Safety Scalpel

- · Safe and multifunctional usage
- Retractable blade for needlestick safety
- before and after usage
- Intuitive pen-like activation
- Protected blade to cut the PEG thread



Safety Puncture Cannula

- Trokar with protective cap increases safety after puncture and during waste disposal
- · Air valve in puncture cannula reduces air escaping from the stomach during endoscopy

Safe ENFit connectors

- Dedicated to enteral application
- Avoids misconnections

Safe and convenient hard blistered tray

- Easy to remove components
- Ready to hand
- Sterile contents fixed by cover during aseptic provision

Safe and convenient documentation

- Patient brochure to educate handling for patients
- 4 sticker with batch number and product name to eases correct documentation
- · Patient implant card which enables and facilitates tracking
- Instructions for use for healthcare professionals



Common pull method for placement

- Easy fixation of PEG tube with thread
- · Radiological control possible
- Easy to combine with an Intestinal Tube for jejunal feeding e. g. gastroduodenal motility problems
- Possible changeover to a balloon system, like Freka Belly Button



Easy aftercare

- · Compact and exchangeable ENFit connectors
- Transparent tube material allows fluid control
- Colour coded accessories make component exchange easier
- Tube clamp prevents reflux of stomach contents

Reliable product quality

for optimal patient comfort, tolerance and safe feeding

- Tissue-friendly polyurethane tube
- · Soft and flexible fixation plates made of silicone rubber
- Latex-free, PVC-free
- Proven for long durability for low exchange frequency*
- Safe enteral connector ENFit, according to the international safety standard ISO 80369-3
- ENFit System provides a locking function to avoid accidental disconnection
- Luer-free IV-incompatible connection System ENFit avoids misconnections

* Siau K, Troth T, Gibson E, Dhanda A, Robinson L, Fisher NC. How long do percutaneous endoscopic gastrostomy feeding tubes last? A retrospective analysis. Postgrad Med J. 2018 Aug;94(1114):469-474. doi: 10.1136/postgradmedj-2018-135754. Epub 2018 Jul 24. PMID: 30042184.



