

## Content of the Freka® Belly GastroTube kit:

- Freka Belly GastroTube, ENFit
- Syringe for balloon filling;  
the syringe size depends on the  
FR-size of the tube
- Instruction for use
- Patient Information
- Patient sticker



## Freka® Belly GastroTube, ENFit

Percutaneous gastrostomy feeding  
tube for long-term enteral nutrition  
in an established stoma tract.

## Ordering Information

FR size	Color identification	Connector	Balloon size	Packaging	Article code
FR 12	White	ENFit	5 ml	1x1	7751704
FR 14	Green	ENFit	5 ml	1x1	7751705
FR 16	Blue	ENFit	10 ml	1x1	7751706
FR 20	Purple	ENFit	15 ml	1x1	7751708

© Fresenius Kabi AG. Reproduction and distribution - only with prior written authorization  
(01.25/AC)



# Freka® Belly GastroTube, ENFit


Percutaneous gastrostomy feeding tube for long-term enteral nutrition in an established stoma tract.



The Freka® Belly GastroTube is designed to support patients in need of long-term enteral feeding or gastric decompression.

Suitable for both hospital and home care settings, it provides a solution for delivering essential nutrition.

Healthcare professionals and caregivers can rely on the Freka® Belly GastroTube which ensures a high standard of care for their patients.



## Who benefits from it?

The Freka® Belly GastroTube is intended for patients who:

- Have an established stoma channel (e.g. after PEG placement) and
- Require long-term enteral feeding directly into the stomach or
- Need decompression or drainage of the stomach



## Who can use it?

- Suitable for both adults and children over 1 year of age
- Available in various FR sizes to meet different patient needs
- Designed for stomata with different tracts lengths



## Which features make it better for patient comfort?

- Made from skin-friendly silicone rubber for gentle, long-term use
- An easy to open external fixation plate keeps the tube securely in place while allowing for comfortable cleaning
- The angled tube position guided by the retention plate makes it more comfortable for the patient
- A built-in tube clamp helps prevent reflux of stomach contents while attaching the Giving Set or syringe.
- X-ray visibility is ensured with a radio-opaque stripe in the tube tip
- FR sizes are easily identifiable by valve color, matching the corresponding color of the Freka® PEG Pro range